
bioquell.com/BI

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Instructions

Bioquell
 Biological
 Indicator

Hydrogen peroxide vapour
 biological indicators

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Bioquell Biological Indicator Instructions

Processing

- Place the Biological Indicators in appropriate locations throughout the target enclosure or in previously validated areas.

- The printed side of the indicator pouch is the inoculated side. This side should always be placed face up although the pouch should be positioned such that all surfaces are exposed to hydrogen peroxide vapour to prevent false positive results.

- Once the cycle is complete, the indicators should be removed from the enclosure and processed/cultured promptly.

Culturing

- Culturing should be conducted using aseptic procedures, preferably in a biological safety cabinet or laminar flow cabinet.

- Open the pouch by peeling back the packaging or alternatively by cutting with sterile scissors.

- Using sterile forceps or other aseptic techniques remove and place the carrier into sterile Broth known to promote sufficient growth of the organism (e.g. Tryptone Soya Broth). Positive and negative controls should be carried out for each cycle.

- Incubate the test and control broth tubes for 7 days at 57.5°C (±4.5°C) and follow the interpretation and disposal guidelines below.

Interpretation

- Tubes should be checked for evidence of growth daily for the full 7-day incubation period.

- Growth of the organism is indicated by turbidity (cloudy appearance) of the broth.

- Turbidity indicates one or more spores survived the bio-decontamination process. No turbidity indicates bio-decontamination process was successful.

Disposal

- Treat expired Biological Indicators and cultured broth as biological waste.

- Expiry dates are marked on each unit. Do not use indicators past their expiry date.

- Incinerate or autoclave expired Biological Indicators and cultured broth at 121°C for at least 30 minutes prior to disposal.

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